



Areas of application for the Aktiv Tramp®

- Improved coordination and stabilisation
- Improved motor function and sense of balance
- Offers unique training possibilities: combinations of training
- Intramuscular and intermuscular coordination in conjunction
- Improved coordination and rhythm
- Multidimensional proprioceptive training
- Low-stress and easy on joints
- Early functional, post-operative use is possible
- Active movement for unsure or older patients
- Supports prevention of osteoporosis
- Stabilising the pelvic floor
- Up to 5Hz mechanically swinging mat

Technical data

- Medical device
- Diameter trampoline = 75 cm
- Dimensions H x W x D = 1,65 x 1,10 x 1,24 m
- Weight approx. 165 kg
- Max. load 150 kg
- Patent protected

Options

Aktiv Tramp® Training poster
Aktiv Tramp® Training poster with stand
Aktiv Tramp® Accessories stand
Aktiv Tramp® Workbook „Therapie und Training“

Do you have any questions? We will be delighted to advise you and provide you an individual and non-binding offer.