

# Pelvic Trainer®



## Areas of application for the Pelvic Trainer®

- Impotence and other sexual dysfunctions
- Aconuresis
- Pelvic instability
- Strength training for pelvic floor musculature
- Prevention of possible future pelvic floor musculature dysfunctions



## Technical data

- Medical device
- Dimensions: H x W x D = max. 117,3 x 50 x 94 cm
- Weight approx. 42 kg
- Max load 185 kg
- Height adjustable seat: 45 – 59 cm
- Sensortube with USB cable and software included
- Angle of shelf for notebook adjustable

## Options

- Pelvic Trainer START UP PAKET (patient flyer, pass for 10 sessions, poster and product film)

**Do you have any questions? We will be delighted to advise you and provide you an individual and non-binding offer.**

**FREI AG** Tel. 07661 / 9336 – 0 | Fax: 07661 / 9336 – 50 | [info@frei-ag.de](mailto:info@frei-ag.de) | [www.frei-ag.de](http://www.frei-ag.de)

**FREI SWISS** Tel. 043-811 47 00 | Fax: 043-811 47 01 | [info@frei-swiss.ch](mailto:info@frei-swiss.ch) | [www.frei-swiss.ch](http://www.frei-swiss.ch)

**FREI Nord** Tel. 040 890 61 63 | Fax 040 890 61 649 | [info@frei-nord.de](mailto:info@frei-nord.de) | [www.frei-nord.de](http://www.frei-nord.de)